

Measuring methods available and examples of their applications

¹³C NMR (carbon nuclear magnetic resonance)

For most organic compounds, information on the number and character of the carbon is essential for confirmation and proof of their structure. As the most common carbon isotope ¹²C does not possess magnetic moment, NMR spectroscopy relies on ¹³C whose natural abundance is only 1.1%. This, combined with four times lower gyromagnetic ratio compared to ¹H, results in much lower measurement sensitivity than in the case ¹H NMR and the measurements are therefore more time-consuming. Routine ¹³C NMR spectra are measured usually with broadband proton decoupling and sensitivity enhancement (3x) by NOE effect. In our CF, a very sensitive cryoprobe (175 MHz at 16.45 T) is available for the measurement of ¹³C -NMR, which significantly reduces the measurement times (tens of minutes).

In analogy to proton NMR, ¹³C NMR allows the identification of carbon atoms in an organic molecule just as proton NMR identifies hydrogen atoms. As such, ¹³C NMR is an important tool in chemical structure elucidation in organic chemistry.

From a ¹H -decoupled carbon (¹³C) spectrum useful information about ¹³C chemical shifts is obtained. Alternatively, carbon multiplicity and ¹H-¹³C coupling constants can be determined from the gated decoupled methodology and quantitative measurements can be made from the inverse gated experiments.

Assignment of the ¹³C spectrum is usually performed with the help of additional NMR methods. As a general strategy:

- Carbon multiplicity is determined from 1D APT or DEPT experiments.
- Correlation with ¹H nuclei via ¹J(CH) is achieved from HSQC or HMQC-type experiments.
- Correlation with ¹H nuclei via ⁿJ(CH) is achieved from HMBC-type experiments.

Examples




Fig. 1. Strychnine nitrate in DMSO-d₆ (approx. 50 mM) ¹H -decoupled ¹³C NMR spectrum with assignment. Spectrometer: AVANCE III HD 700, Probehead: 5 mm CPTXO 13C/15N-1H/D with z gradients, Experiment time: 10 min




Fig. 2. Cefotaxime Acid (cephalosporin antibiotic), ¹H-decoupled ¹³C NMR spectrum with assignment, Spectrometer: AVANCE III HD 700, Probehead: 5 mm CPTXO 13C/15N-1H/D with z gradients, Experiment time: 10 min

For quantitative applications, it is necessary to obtain reliable integration values. In this case, the ¹³C spectrum is recorded with a 90° pulse and a long relaxation period (5 * T₁(¹³C)) is used to avoid partial signal saturation. It can be used for example for characterization of synthetic polymers.